

**June 1972**



# **orthotics and prosthetics**





An Official Journal of



**EDITOR**  
Alan A. Smith, Jr.

**TECHNICAL EDITOR**  
A. Bennett Wilson, Jr.

**EDITORIAL BOARD**  
Bert R. Titus, C.P.O.,  
Chairman  
Mary S. Dorsch, C.P.O.  
William M. Brady, C.P.  
Colin McLaurin, Sc.D.  
Robert E. Tooms, M.D.

**Orthotics and Prosthetics:** The Orthopedic and Prosthetic Appliance Journal is issued in March, June, September and December. Subscription price, payable in advance, is \$6.50 a year in the Western Hemisphere. Rate elsewhere is \$7.50 a year. Publication does not constitute official endorsement of opinions presented in articles. The Journal is the official organ of the publisher, The American Orthotic and Prosthetic Association. All correspondence should be addressed to: Editor: **Orthotics and Prosthetics**, 1440 N St., N.W., Washington, D.C. 20005. Telephone, Area Code 202, 234-8400.

# orthotics and prosthetics

THE JOURNAL OF THE  
ORTHOTIC AND PROSTHETIC PROFESSION

VOLUME 26 • NUMBER 2 • JUNE 1972

## Contents

- 1 Editorial—What Is the Role of the  
"Restorations Prosthetist" in Rehabilitation?  
*Maurice A. LeBlanc, C.P.*
- 3 Prosthetic Considerations for the Female  
*Mary S. Dorsch, C.P.O.*
- 6 Experience With the Stuttle Spinal Brace  
*Paul Campbell, M.D.*
- 12 Fabrication of the Stuttle Spinal Brace  
*Dale O. Butler*
- 17 Polyethylene Girdle for the  
Milwaukee Brace  
*Jerome E. Skahan, C.O.*  
*John J. Skahan, C.O.*  
*Walter E. Dase, Technician*
- 23 Use of A Cast-Brace Application Device,  
Plaster Brim and Ankle Support for  
Fracture Bracing  
*Max Lerman, C.O.*
- 31 Use of Photography on Cast Preparation  
*James H. Tyo*  
*George Hall, Jr., C.P.*
- 34 Population and Access to Prosthetists  
*G. E. (Ned) Sharples, Ph.D.*

XXVII Classified Advertisements

# Cut it. Mold it. Bend it. Shape it.

## No other splint material works so easily.

Wrap around any contour. Lend support to any injured body part. With attractive ORTHOPLAST® Splints, the most versatile splinting material available.

- ☐ A lightweight low-heat thermoplastic featuring maximum ability to be molded.
- ☐ Strong, cohesive, can be made to adhere to itself without glue, simply by heating.
- ☐ Can be riveted, strapped, braced, hinged, bonded, butt-bonded to suit all your splint-engineering needs.

(Can also be reheated and remolded to refit an immobilizing patient whose splint needs have changed.)

Educational literature, medical reprints, splint patterns and samples may be obtained by writing to Department J-235, Johnson & Johnson, New Brunswick, New Jersey 08903.

**ORTHOPLAST** Splints

TRADEMARK



*Johnson & Johnson*



# THE AMERICAN ORTHOTIC AND PROSTHETIC ASSOCIATION

## OFFICERS

President—**Mary S. Dorsch**, C.P.O.  
New York, New York

President-Elect—**Claude J. Lambert**, C.P.O.  
New Orleans, Louisiana

Secretary-Treasurer—**Durward R. Coon**, C.P.O.  
Detroit, Michigan

Vice-President—**Robert V. Bush**, C.P.O.  
Albuquerque, New Mexico

Immediate Past President—**Roy Snelson**, C.P.O.  
Downey, California

## REGIONAL DIRECTORS

Region I—**Siegfried W. Paul**, C.P.O.  
Newington, Connecticut

Region VII—**Dale E. Clark**, C.P.  
Waterloo, Iowa

Region II—**JACK GOLD**, C.P.O.  
Orange, New Jersey

Region VIII—**JERRELL BALLARD**, C.P.,  
Amarillo, Texas

Region III—**Robert E. Teufel**, C.P.O.  
Elizabethtown, Pennsylvania

Region IX—**Herman Hittenberger**, C.P.O.  
San Francisco, California

Region IV—**HENRY G. BOWDEN**, C.O.,  
Warm Springs, Ga.

Region X—(Acting) **Leonard W. Jay**, C.P.  
Salt Lake City, Utah

Region V—**Wade L. Barghausen**, C.P.  
Columbus, Ohio

Region XI—**William E. Brownfield**, C.P.  
Boise, Idaho

Region VI—**Loren D. Jouett**, C.P.O.  
Chicago, Illinois

## AMERICAN BOARD FOR CERTIFICATION IN ORTHOTICS AND PROSTHETICS, INC.

President—**J. D. Ferguson**, C.P.O.  
Chapel Hill, North Carolina

President-Elect—**Donald F. Colwell**, C.P.  
Grenada Hills, California

Immediate Past President—  
**Robert E. Fannin**, C.O.  
Columbus, Ohio

Vice President—**Alan G. Smith**, M.D.  
Brooklyn, New York

**Robert C. Hamilton**, M.D.  
Chicago, Illinois

Secretary-Treasurer—**Robert F. Hayes**, C.P.  
Pittsfield, Massachusetts

**Robert E. Tooms**, M.D.  
Memphis, Tennessee

**Fred J. Eschen**, C.P.O.  
New York, New York

## AMERICAN ACADEMY OF ORTHOTISTS AND PROSTHETISTS

President—**William M. Brady**, C.P.  
Kansas City, Missouri

President-Elect—**Arthur W. Guilford, Jr.**, C.O.  
Downey, California

Vice-President—**Henry F. Gardner**, C.P.O.  
New York City

Secretary-Treasurer—**Daniel G. Rowe**, C.P.O.  
St. Paul, Minnesota

Immediate Past President—**Ralph R. Snell**, C.P.  
Memphis, Tennessee

### DIRECTORS

**H. Blair Hanger**, C.P.  
Park Ridge, Illinois

**Hugh J. Pantan**, C.P.  
Miami, Florida

**Charles L. Jones**, C.P.  
Hambra Heights, California

**Herbert B. Warburton**, Executive Director  
**Alan A. Smith, Jr.**, Asst. Executive Director

## Index To Advertisers

**JUNE 1972**

<b>C. H. ALDEN SHOE CO.</b>	Inside Back Cover
<b>BECKER ORTHOPEDIC APPLIANCE CO.</b>	XXII
<b>S. H. CAMP AND CO.</b>	VIII
<b>C. D. DENISON CO.</b>	XXIII
<b>FILLAUER SURGICAL SUPPLIES, INC.</b>	XVII
<b>Fiway Manufacturing Co.</b>	XVI
<b>FLORIDA BRACE CORP.</b>	XXV
<b>FREEMAN MANUFACTURING CO.</b>	XIV
<b>Herbst Shoe Mfg. Co.</b>	V
<b>Johnson &amp; Johnson</b>	II
<b>James R. Kendrick Co.</b>	XV
<b>KINGSLEY MANUFACTURING CO.</b>	X
<b>KNIT-RITE, INC.</b>	XVIII

<b>L. LAUFER &amp; CO.</b>	XXIV
<b>Major Lab Manufacturers</b>	XXIV
<b>Orthopedic Splints, Inc.</b>	XI
<b>PEL SUPPLY CO.</b>	XII
<b>POLYCADENCE, INC.</b>	VI, VII
<b>REALASTIC INDUSTRIES</b>	XXII
<b>Robert O. Porzelt Co.</b>	XXVI
<b>Roden Leather Co.</b>	IV
<b>E. J. SABEL CO.</b>	XIII
<b>TRUFORM ANATOMICAL SUPPORTS, INC.</b>	XX, XXI
<b>WILH. JUL. TEUFEL</b>	XIX
<b>WASHINGTON PROSTHETIC SUPPLIES</b>	IX

Advertisers shown in bold-face type are members of The American Orthotic & Prosthetic Association.



**SERVICE  
OUR  
BYWORD**

—Exclusively—

**LEATHERS FOR THE  
ORTHOPEDIC and PROSTHETIC TRADE**

**LIMB & BRACE LEATHERS:**

- Orthopedic Horsehide
- Orthopedic Cowhide
- Orthopedic Elkhide
- Glazed & Molding Cowhide
- Calf—Kip—Sheep—Pigskin

**HOSPITAL LEATHERS:**

- Shearlings for Pads
- Carving & Tooling
- Molding
- Russet Strap

*Satisfied Customers Throughout the United States Since 1924*


**ORDERS SHIPPED  
THE SAME DAY**

**RODEN LEATHER CO., INC.**  
1725 CROOKS ROAD  
ROYAL OAK, MICH. 48068  
Area Code 313-542-7064



# Control can be more comfortable

with the new  
**Child Life**  
pre-walker  
abduction  
boot



The Sof-Top seamless back and the serrated tongue assure double protection against chafing. These features are available on Child Life abduction and straight last boots... and there is *no charge for our Mismatch Service*. All Child Life Specialty Boots are available in sizes from 000 to 4, in narrow and medium widths.


**Two-Strap Abductor Boots**  
**Single-Strap Abductor Boots**  
**Advanced Straight Last Boots**  
**Basic Straight Last Boots**

- 1 Sof-Top Seamless Back...** prevents chafing and blistering while retaining the high back necessary to hold the heel securely in the shoe
- 2 Serrated Tongue...** prevents irritation to the instep area
- 3 Wider Toe Opening...** gives you a better view of foot position within the shoe



Herbst Shoe Manufacturing Company  
2775 N. 32nd St., Milwaukee, Wisconsin 53210


*This amputee has a hip disarticulation equipped with an HD-11 Polycadence hip joint and an AK-4 knee with a SACH foot\**



## **Restore confidence to your seriously handicapped patient**

The Polycadence AK-4 has given new hope to many. Here is a proven, case-tested knee that gives confidence to the most cautious patient. The polycentric AK-4 joint has built-in control at every position. This firm support is instantly felt by the patient. He can sense the limb's response. He can hear the slight sounds of position change. And he feels the ease of clearing the foot at each step...without undue elevation. This gives the patient an early success in balancing and a definite psychological lift. From the moment of fitting these qualities are discernible—in contrast to the ordinary over-cautious steps of the seriously handicapped amputee.





*During extension the built-in alignment of forces locks the knee joint into a firm position without patient adding pressure.*

*During flexion the polycentric balance automatically shifts to allow joint movement without fear of toe drag.*

The Polycadence AK-4 is mutually beneficial to both you and your patient:

**Alignment** Simple to work. It's built with a 2° forward tilt to the pylon, pre-aligned, and requires no further refined adjustments.

**Versatile** The Polycadence modular system is uniquely adaptable. The same unit is used for temporary fitting and later for permanent prosthesis.

**Durability** Through years of testing, a unique selection of materials has been chosen for trouble-free, long-lasting performance. Over 500 units presently in use with the seriously handicapped attest to this.

**Stability** The knee's built-in controllable support provides the patient with immediate confidence of his stability.

**Flexibility** Complete ease and control of movement from one position to another. Foot clearance without undue elevation. (See illustrations at left.)

**Cosmesis** Pre-shaped covers for men. Shapeable foam block for women for custom matching of existing limb. Caucasian or Negroid.

Many types of patients can benefit from the AK-4: Bilateral involvements, short stumps, hip disarticulations, hemipelvectomies, geriatric, the mentally impaired... It is no longer a question of whether they are able to walk. It becomes, instead, a challenge to both you and your patient as to *how many* accomplishments are possible for him. No longer fearing for his balance, he can concentrate on mastering new achievements.

# polycadence AK-4 PATENT PENDING


POLYCADENCE, INC., 6070 WEST 55th PLACE,  
ARVADA, COLORADO 80002, PHONE 303/423-3412.

**POLYCADENCE AK-4 IS A V.A. CONTRACT ITEM.**

For complete technical information and price list, write or call the above. All literature and sales of materials are available to Certified Prosthetic Facilities only. No exceptions.

*\*The HD-11 Polycadence hip joint is another progressive development of the Research Department of Polycadence, Inc., and further information is available upon request.*

One of these elastic stockings provides maximum compression at the ankle, gradiently decreasing to calf and thigh—to provide proper support without downward compression.



And was fitted from stock.

The basic problem with most elastic stockings has been that they don't provide the required compression characteristics at the ankle and throughout the leg. Or, they are capable of compression but must be hand-made to provide it. Or, they provide compression initially but lose it within a few short weeks. Or, they are so cosmetically distasteful the patient won't wear them.

None of these problems exist in the clinically approved Sigvaris® stocking now being distributed in the United States by S. H. Camp and Company.

Available in calf, half-thigh, and half-thigh plus lengths; in two compression strengths, normal (30-40 mm Hg) and super (40-50 mm Hg) and in small, medium and large sizes. Seamless, easily applied, and scarcely noticeable when worn under regular hose, Sigvaris stockings are readily accepted by patients. A basic stock enables you to properly fit over 85% of all patients — immediately.

® Registered Trademark of Ganzoni & Cie Ltd., St. Gall-W, Switzerland



Distributed by

**CAMP**





# HAVING PROSTHETIC PRODUCTION PROBLEMS?

TOTAL CENTRAL FABRICATION IS NOW AVAILABLE AT OUR NEW MODERN PROSTHETIC LABORATORY. YOU CAN QUICKLY OBTAIN THE HIGHEST QUALITY UPPER AND LOWER EXTREMITY APPLIANCES READY FOR TRAIL FITTING AS WELL AS COMPLETE FINISHING SERVICE. MANY PROSTHETISTS HAVE FOUND REAL BENEFITS IN CONVERTING TO CENTRAL FABRICATION — HOW ABOUT YOU?

WE ALSO CARRY A FULL LINE OF PROSTHETIC COMPONENTS AND SUPPLIES READY FOR IMMEDIATE DELIVERY. OUR INVENTORY INCLUDES THE "MARVEL" LINE OF KNEE SHIN ASSEMBLIES AND SACH FEET, FOR WHICH WE ARE AN EXCLUSIVE DISTRIBUTOR. COMPARE OUR PRICES . . . . .

SEND FOR YOUR FREE CATALOG TODAY . . . . .



## WASHINGTON PROSTHETIC SUPPLIES

40 PATTERSON STREET, N.E.  
WASHINGTON, D.C. 20002

Telephone (202) 628-1037

• **PRESENTING** •

**The KINGSLEY High Profile Molded SACH Foot**

**Approved by the Veterans Administration**

**Shaped to the VA Models**



Consistent toe resistance.

Triple density heel.

Can be shaped to half sizes and narrow widths.

Natural flesh color.

Completely odorless.

Available in extra soft, soft, regular, firm and extra firm heel cushions  
(toe stiffness remains the same!).


Foam extends over top of keel.

Stress relief at the top of the bolt hole.

**Kingsley mfg. co.**

1984 PLACENTIA AVENUE • COSTA MESA, CALIFORNIA 92627 • PHONE 714-548-6116






An advanced design of the Denis Brown night splint. Anodized aluminum bars in sizes from four inches to thirty-six inches with your logo embossed on the bar if desired. Now offered with skate-clamp attachment, incorporating a full protractor. Plates are available in two sizes. For same day shipping contact

ORTHOPEDIC SPLINTS, INC. — P.O. BOX 508  
LYNBROOK, NEW YORK 11563

Or phone (516) 887-9280

. . . DAY AFTER  
DAY AFTER DAY.





**PROSTHETISTS:**  
**AMP-AID can be of  
great help to you  
and the patient —  
also the therapist!**

Gait training of a new amputee is often delayed or slowed down due to a red or discolored area on the stump. The cause is usually that the stump is a little swollen that day or the wrong thickness of stump sock is being used. However, the amputee or therapist doesn't know this or they aren't sure, so training is stopped or delayed. If it happens again, the patient is referred back to the Prosthetist for a check up and possible adjustments to the prosthesis.

AMP-AID ointment may be a solution to this problem. Apply about 1/16" thick over the red area and roll the stump sock on. The red area slides in the AMP-AID, which is not absorbed into the skin, but stays on the surface and acts as a lubricant. Training can be resumed, with caution, so you and the therapist save time and effort—the amputee saves time, money and avoids discomfort. Make sense? Worth trying? We think so!

You'll also find AMP-AID a most dependable ointment for treating many other skin irritations such as minor burns, chafing from orthoses, shoes, tops of ski boots, etc. It is safe to use, with no side effects. Registered with the Federal Food and Drug Administration.

AMP-AID contains 95% active ingredients: Olive Oil, Cottonseed Oil, Stearine, Petrolatum, Camphor, Phenol, Boracic Acid, Zinc Oxide and Lanolin.

AMP-AID may be obtained from your local Prosthetist or Orthotist.

## **PEL SUPPLY COMPANY**

4472 WEST 160TH ST., CLEVELAND, OHIO 44135 — PHONE: 216-267-5775

PAUL E. LEIMKUEHLER, PRESIDENT

Guess what Sabel-Stride Rite did to  
eliminate rubbing and chafing  
at the heel?

They cut out the problem!

Sabel-Stride Rite  
did it by removing  
the non-essential  
leather at the top  
of the heel that  
causes irritation when  
unnatural pressure is  
applied to the foot  
during splinting  
or after casting.  
Sabel found the way.  
With a little less shoe,  
& a lot more comfort.



**SABEL Division**  
The Green Shoe  
Manufacturing Co.  
makers of  
Stride Rite Shoes

For information, write:  
E.J. Sabel & Co.,  
Benson-East, P.O. Box 644,  
Jenkintown, Pa. 19046





# Freeman dealers don't need a warehouse...



they have ours at their finger tip

In the United States there are four strategically placed Freeman warehouses. All are fully stocked with the Freeman line of high quality surgical support garments. By calling the nearest warehouse collect, any Freeman dealer can have his order processed and shipped within 24 hours, or less. Thus Freeman dealers need only basic inventory...freight costs are less, too. *Write for complete information and the new Freeman catalog.*

#### **FREEMAN WAREHOUSES:**

Ephrata, Pennsylvania 17522, P.O. Box 277, phone 717-733-4261

Sturgis, Michigan 49091, P.O. Box J, phone 616-651-2371

Anniston, Alabama 36201, P.O. Box 1791, phone 205-237-0611

West Covina, California 91790, 636 S. Sunset St., phone 213-338-1618

*freeman*

FREEMAN MFG. CO.

Box J, Sturgis, Mich. 49091



# DEMAND THE FINEST FOR YOUR CUSTOMERS

## KENDRICK'S dorso-lumbar supports



Kendrick Dorso Lumbar Support No. 2030 features 21" back depth, an elastic release at the groin and a hook and eye opening at the left side.

### **PRESCRIBED BY PHYSICIANS...SOLD WITH CONFIDENCE**

For your customers there should be no substitute for the finest. Kendrick's Dorso Lumbar Support for women offers extremely beneficial protection for the dorsal region as well as the lumbar spine. The Dorso Lumbar Support can be worn following fractures, in cases of faulty posture defects such as scoliosis, and for arthritis of the spine.

Model No. 2030, illustrated, is available in stock sizes from 30 to 44, hip measurement — hip development, 6"—8". The support is well boned in the back with two 18" shaped steels at the center to brace and support the entire back. Three traction straps permit variable pressure. Support through the dorsal region is increased by the shoulder straps which are looped under the arm and joined in the front. Also available in Dacron Mesh Style No. 8430.

Specify surgically correct Kendrick products for your customers. Their confidence in you will result in more frequent repeat sales and higher profits.

*You can recommend Kendrick with Confidence*

**JAMES R. KENDRICK COMPANY, INC.**  
Philadelphia, Pa. 19144      New York, N. Y. 10016

***Kendrick***  
S I N C E 1 8 5 3

# TAYLOR - KNIGHT - WILCOX

## OR COMBINATION OF WILCOX-TAYLOR - OR - KNIGHT-TAYLOR

### SINGLE - OR - DOUBLE PELVIC BAND STYLE



**TAYLOR**  
straight pelvic

All parts are interchangeable  
Use any combination of bands

3 style single pelvic - 2 style thoracic bands  
Plus

The hard to make double pelvic  
Add cervical support—chest plate—auxiliary straps  
Assemble any style in matter of minutes  
In your own shop - no special tools required



**TAYLOR**  
double pelvic



**KNIGHT**  
Downflare  
Thoracic  
Butterfly Pelvic



**KNIGHT-  
TAYLOR**  
Straight  
Thoracic  
Upflare Pelvic



The  
Modern  
Way  
To  
Assemble  
Spinal Braces  
All  
Made



**KNIGHT-  
TAYLOR**  
Double-Pelvic  
Cervical  
Support



**KNIGHT**  
Straight  
Thoracic  
Double Pelvic

Possible with the use of  
**PRE-FABRICATED PARTS**

All fitters must be instructed on the assembling procedure of these braces  
There is no other way this product can be purchased  
Enroll now in a class closest to your area  
1972

#### LOCATIONS OF CLASSES FOR THE NEXT 3 MONTHS

**SEE  
HOW  
OUR  
PRODUCT  
WORKS  
YOU BE  
THE JUDGE**

June	July	August	September
3rd Chicago, Ill.	8th Hartford, Conn.	5th San Fran- cisco, Ca.	9th Denver, Colo.
17th Cincinnati, Ohio	22nd Milwaukee, Wisc.	19th Seattle, Wash.	23rd Minneapolis, Minn.
—	—	26th Pittsburg, Pa.	—

**WRITE  
FOR  
COMPLETE  
SCHEDULE  
OF  
SEMINARS  
AND  
APPLICATION**

516 S. Howard Ave.

**FIWAY MANUFACTURING CO.**

Tampa, Florida 33606



# What is the Role of the "Restorations Prosthetist" In Rehabilitation?

Appearance is important to most people—men and children as well as women. It is infrequently discussed in clinics but is frequently in the minds of patients.

People in need of cosmetic restoration include amputees, various orthotic patients, women with mastectomies, those with facial disfigurements, and many others. Whether it is minor such as a missing finger or major such as a missing arm, successful cosmetic treatment enhances body image and makes people happy.

To provide acceptable cosmesis is a difficult job, for it includes not only static appearance but how something looks in motion and how it feels, sounds, and smells. An artificial eye is much more effective if it moves with the good eye. An upper-limb prosthesis appears much more natural if it swings, feels soft, and makes no mechanical clicking noises.

By and large, most of the prosthetic and orthotic facilities are more occupied and better trained to provide functional restoration

rather than cosmetic restoration to patients with artificial limbs and braces. There are several maxillofacial centers whose work involves facial buildups and replacements of the nose and ears; most are related to dental schools and clinics. Twelve Veterans Administrations hospitals and the VA Prosthetic Center provide cosmetic care for VA patients. And finally, there are a small number of individuals practicing privately to provide cosmetic service.

One manufacturer supplies most of the gloves, leg covers, and maxillofacial materials being used. This manufacturer has repeatedly tried to avoid dealing directly with patients. The company has given courses, written instructions, and provided kits to allow cosmetic work to be done in the field. However, so many people in the field are so busy providing function and just trying to stay abreast of new developments that they neither have the time nor desire to do the rather meticulous and time consuming work of providing cosmesis. The result is that many patients end up knocking at



the manufacturer's door or being shuffled to the few people around who perform this service. In general, the present situation seems to be that patients are obtaining cosmetic treatment one way or another, but that more can be done for more people.

The "restorations prosthetist," or "anaplastologist" as coined by Mr. Peyton Massey, is an individual who provides cosmetic restoration service for patients. There are only a handful of them yet there are many patients who can benefit from this "added" service. In many respects, he is to rehabilitation as the industrial designer is to household products; he beautifies the product and enhances its use for the consumer.

Why not put the restorations prosthetist in the rehabilitation center so his services are available to patients as part of a comprehensive treatment program? He would certainly get referrals from prosthetic, orthotic, dental, and eye clinics and from various surgeons. It probably would not be a large service in numbers but could be a very important one.

There was a Negro man who refused his below-elbow prosthesis

three times under the guise of mechanical difficulty until someone had the idea of making it a shade lighter in color. There are several young post-polio girls with atrophied legs who walk more confidently now with leg buildups. There are people with prosthetic eyes, noses, and ears who can at least face the public even though function has not been restored. These things are what the restorations prosthetist can do.

If the discipline of "anaplastology" were added to the rehabilitation center, it would help to bring together the scattered efforts in cosmetic restoration; it would relieve many prosthetists, orthotists, and dentists of work they would just as well not bother with; and it would give a big psychological boost to patients in serious need of it. And if it were really successful, a new allied-health specialty could be developed with its own education, training, and qualification program. It seems like a possibility at least worth exploring.

MAURICE A. LEBLANC, C.P.

Committee on Prosthetics  
Research and Development  
National Research Council

# Prosthetics Considerations for the Female<sup>1</sup>

Mary S. Dorsch, C.P.O.<sup>2</sup>

I am very happy to have been invited to participate in this symposium on patient management. I do feel that this is an area that has been completely neglected in our professional thinking. In striving to perfect the technology, we have forgotten the patient—and above all—the FEMALE patient. It has been said that there is no difference between the male and female patient. Well, to that, I say 'Vive Le Difference', and will discuss now the problems of the female patient.

To all men, women are an enigma, but no greater problem has any man, than one who must try to rationalize amputation to a female

patient. Understanding of the emotional and physical problems of the female patient can influence her adjustment and acceptance of the prosthesis. In many instances, we may eliminate some of the problems by approaching the situation with a little more finesse and better understanding.

One of the very important things that we must not do—that is not to de-personalize the patient. Do not label her—do not call her 'an above-knee amputee' or 'a below-knee amputee' or a 'hip-disarticulation' case. There is no need for that because the patient is a whole being, and it is very important to consider this. Total rehabilitation requires planning by the entire team, beginning with the surgeon and terminating with the fitting of the prosthesis.

The following considerations will be a great help to the patient and to the members of the team. First, a

<sup>1</sup> This paper was presented at the Amputee Management Symposium sponsored by the Amputee Center, Rancho Los Amigos Hospital, University of Southern California, in association with the American Academy of Orthotists and Prosthetists and the National Academy of Sciences, March 20-22, 1972.

<sup>2</sup> 109 East 29th Street, New York, N.Y. 10016.



woman is a woman, then a patient. Any indication that she is considered only an 'amputee' would be an emotional and psychological setback to her vanity. The physician, in trying to alleviate the initial shock of the amputation, will describe the subsequent prosthesis in such a manner that the prosthesis cannot meet the mental image in the patient's mind. Although the modern prosthesis is an excellent substitute, functionally and cosmetically, for the missing member, there are limitations within which the prosthetist must work.

In my conversations with physicians, many have told me of their unawareness of new developments in our profession. Therefore, I firmly believe that the prescription for a prosthesis should be a cooperative venture between the physician and the prosthetist. The physician being responsible for the medical aspects, the fitting or fashioning of the prosthesis is the responsibility of the prosthetist, and that is where his or her expertise will prevail.

Psychological setback is seen in some surgeons as well as in the patients. And it may surprise you to learn that some surgeons equate amputations to death. (They are not alone in this feeling as insurance companies equate double amputations to death.) However, with immediate post operative fitting procedures, as well as early fittings, I do believe a closer follow-up by the surgeon will relieve him of this frustration and help him to realize he has given the patient the ability to return to being a normal human and leading a useful life.

Many of us in the profession rarely stop to consider the mental

image our words convey. Try not to use the word 'stump' which is a degrading and painful word. It has been said that a picture is worth a thousand words, and yet, there are just two words that conjure up most gruesome pictures to persons who have lost limbs. It is our responsibility to replace these words, as well as the limbs. The words are 'amputee' and 'stump'. While these words affect men severely (although they are supposed to be the hardier of the species, in spite of the Women's Liberation Movement), the loss of a limb becomes a traumatic experience to a woman, and not only for vanity's sake.

Little girls are likely to have fewer anxieties and vanities, so they are just great to work with. However, when you go into the teen-age stage, you really have a problem in their manner of dress and how their friends, male and female, are going to accept them.

Then from the teen-ager, we go into the adult stage and you have a young woman who is planning courtship, marriage, and motherhood. There are many problems to be considered at this stage. I wonder how many doctors or prosthetists here have had a young girl request your advice regarding pregnancy. The first trimester does not present too many problems, but during the second trimester, adjustments for comfort and control are necessary. I feel that we should be more aware of these very personal situations, and, of course, we should have the answers.

The below-knee prosthesis presents little problems. However, prostheses for higher amputations and



hip disarticulation cases require adjustments to suspension and refashioning of the socket in the anterior aspect, with perhaps a leather replacement from time to time until confinement.

At the time of confinement, the prosthesis should remain on the patient, as this will permit better positioning on the delivery table and also comfort the patient both physically and mentally.

Now, as you consider higher levels of amputations, this will have to be evaluated by the obstetrician.

We are in a great era regarding clothes for the females—the skirts are long—pant-suits are in—and the shoes are great. However, speaking of shoes, do fit the young ladies, and the not-so-young ladies, with several feet at the initial fitting—that is, high heel, low and flat heels. This will eliminate a bilateral frustration. They will have a selection of shoes and you will have a happier patient, one who is satisfied because she realizes that you fully understand her situation.

I respect whatever procedures you prefer—that is the immediate post-surgical procedure—or a rigid dressing with early fitting—as long as the patient leaves the hospital wearing a preparatory prosthesis with good cosmetic effect. This is true for lower

limb cases as well as for upper limb cases. We have been fitting modular prostheses for upper extremities for the last thirty years. The interscapulothoracic amputation is very disfiguring, and a prosthesis should be fitted before the patient goes home.

The measurement and cast can be made prior to the amputation, and final fitting can be carried out when the sutures are removed.

We advise that you fit the patient with an interscapulothoracic joint as well as a complete prosthesis. When relaxing at home this arrangement will be a little more comfortable and will give a better cosmetic effect since it supports clothing in a presentable manner.

I sincerely hope what I have said will answer some little problems you may have. I also hope it will arouse some questions.

In closing, keep this in mind, please—limbs suggest parts of beautiful trees as well as bodies, and in a measure block out the horrendous bloody picture of a sawed off part of one's body. Trees lose limbs through storms or judicious pruning, yet we never refer to a tree as having an amputation. Treat patients as gently with words as you do with your actual treatment.

Try hard—Thank you.

# Experience With the Stuttle Spinal Brace

Paul Campbell, M.D.<sup>1</sup>

The problem of lumbosacral disease is universal after the age of twenty, and a fairly high percentage of those involved at one time or another require orthopedic management. Most of these patients can be controlled without bracing, but a significant percentage have enough instability or degenerative changes in the "low back" to require the use of external support.

The history of back bracing goes back many centuries but became somewhat scientific approximately only about a hundred years ago with the development of the Taylor "Spinal Assistant" (Fig. 1 & 2). The Taylor brace is effective and is still used commonly. Various modifications have been made including the Knight and Goldthwait braces. A shorter version that extends only to the lower thoracic region is the most

commonly used modification, called either a chair-back brace or cage brace. Another commonly used brace that produces hyperextension is the Jewett brace. All of these devices are rigid in construction and attempt to immobilize the lumbosacral region by grasping the pelvis below and a portion or all of the thoracic regions above by encirclement of the chest. All of these braces are ineffective to some degree because the pelvis cannot be grasped adequately to prevent motion in the sagittal plane of the two lower joints in the lumbosacral spine in movements like sitting and bending forward from a standing position.

A definitive study of these types of braces was made by Norton and Brown in the early fifties (1). A few of their comments are worth repeating because they elucidate this problem. They stated that the effective-

---

<sup>1</sup> Portland Orthopedic Clinic, Portland, Oregon.






Figure 1  
The Taylor "Spinal Assistant" (from 2).

ness of a brace is related to the patient's comfort in using a brace. In other words, a brace that is uncomfortable in use will not be effective because it will not be worn snugly enough to be effective. They further commented that "if the object of treatment is to immobilize the lumbosacral spine, either sitting must be avoided or the apparatus must be effective in the sitting position." Since our culture involves sitting in most of our daily activities including eating, transportation, working and recreation, it is obvious that an effective brace must function in the sitting position. It is in this position that the rigid braces mentioned before are ineffective because they are unable to adjust to the change in relationship between the pelvis and spine when the position of sitting is changed to standing. The third comment that is pertinent is, "It seems highly unlikely that any device ap-

plied to the exterior of the body can effectively splint the lumbosacral region." This is the significant problem facing those prescribing or making "low back" braces.

In an attempt to solve this problem Williams developed a functional brace incorporating a hinged bar as the fixation device across the lower thoracic spine region and a fixed bar beneath the buttocks with the third point of fixation being a broad apron across the abdominal region. (Fig. 1 & 2). The hinge in this brace allowed it to accommodate to the changes in posture from the standing to the sitting position, maintaining the lumbosacral and the L4-5 joints in a flexed position throughout. The Williams brace incorporates lateral uprights to prevent rotation and lateral bending, thus affording the maximum degree of rigidity while still allowing essential function.

Experience with the Williams







Figure 2  
The Williams Brace (from 2).

brace was disappointing in that the patients complained of pressure in the lower thoracic spine region and tended to either wear the brace too loosely or to discontinue its use against advice. Occasionally a patient would return with the brace applied in an upside down position stating that despite the fact the brace did not fit too well in such a position, it was more comfortable and could be tolerated better in an upside down position. Our practice was to admonish the patient to wear the brace properly and to follow the instructions of the surgeon and orthotist.

Dr. Fred Stuttle of Peoria, Illinois, evidently had similar experience with patients, but was wise enough to listen to the patients and to respond to the strong hint that perhaps improvements in the design of this brace were possible. It was a very simple modification to

change the arrangement of the posterior bars so that the hinged bar was placed across the well-padded gluteal region and the fixed bar across the less protected lower thoracic spine area. No other significant modifications were made to the Williams brace by Stuttle who presented his version (Fig. 3) in a scientific exhibit at the 1962 Annual Meeting of the American Academy of Orthopaedic Surgeons. Because of the obvious advantage in function of the Stuttle brace compared to the Williams design, we adopted this brace for use in our clinic.

For the last ten years we have used more than 700 Stuttle braces and are continuing to use approximately three each week for treatment of low back disease. Minor modifications in design have been made by our chief orthotist, Mr. Dale Butler, but the essential functional ele-




Figure 3  
The Stuttle Brace.



ments as first devised by Dr. Williams and modified by Dr. Stuttle have not been changed. Patients almost uniformly tolerate this brace well and continue to use it as long as prescribed. In fact, occasional patients have continued to use the brace even after they were advised to discontinue it because they felt more secure with the brace and could see no reason to discontinue it. There have been no serious complications in using the brace and only occasionally have patients been bothered with irritation of the lower lumbar spine requiring a more rigid type of fixation. This brace is not effective in treating the upper lumbar spine area; it is only effective for the lower two joints of the lumbosacral spine. These two joints are normally hyperextended owing to the pelvic tilt, and thus are best protected by bracing in flexion, because straightening the lordotic curve lessens the sheer component of the forces applied.

We do not consider this brace as corrective of postural abnormalities even though it has been recommended for such use. This brace functions by immobilizing the spine in a flexed position, but does not carry a spine into a greater range of flexion than was possible prior to the application of the brace. The immobilization in a flexed position has the added advantage of tending to open the intervertebral foramina when used for immobilization to achieve spinal fusion, thus affording maximum room for the nerve roots to pass from the spine. The fact that this brace acts to maintain flexion of the lumbosacral spine in all body positions helps it achieve maximum

immobilization as the patient's posture changes from standing to sitting and vice versa. It is obvious that the less sitting a patient does while wearing the brace, the less motion will occur at the involved joints and so our instructions to patients in whom we hope to achieve spinal fusion is "avoid sitting as much as possible during the immediate post-operative period."

## CONCLUSION

We have such confidence in the efficacy of this type of low back brace that we use it for many types of diseases of the lumbosacral spine region where immobilization in a functional position is indicated. Where less stringent support is required we still use a canvas type of low-back support, usually incorporating a pad to fit into the lumbar lordosis. Occasionally for compression fractures in this region we use a hyperextension brace or a Jewett type of back brace. For all other conditions requiring external immobilization, the Stuttle brace is our choice. Specifically we use this brace after spinal fusions, after lumbosacral laminectomies, or for those patients treated conservatively who require more rigid support than the reinforced canvas support can afford.

It can be concluded that a brace to immobilize the lumbosacral spine will be effective only if it is firmly applied, only if it is used continuously, and only if it is tolerated long enough to be effective. We believe that the Stuttle modification of the Williams flexion brace makes it possible for these requirements to be fulfilled.



## REFERENCES

1. Norton, Paul L., and Thornton Brown, The immobilizing efficiency of back braces, *J. Bone and Joint Surg.*, 39A:1:111-138, Jan. 1957.
2. Thomas, Atha, Appliances for the Spine and Trunk, pp. 179-250, *Orthopaedic Appliances Atlas*, Vol. 1, J. W. Edwards, Ann Arbor, Mich. 1952.
3. Williams, P.C., Lesions of the lumbosacral spine, part II, *J. Bone and Joint Surg.*, 19:690, 1937.

# Fabrication of the Stuttle Spinal Orthosis

Dale O. Butler<sup>1</sup>

In 1962, Dr. Fred Stuttle presented to the American Academy of Orthopaedic Surgeons a new concept of flexion-type spinal orthoses, utilizing the true Williams pattern of back bracing, but reversing the thoracic and pelvic band and placing the pelvic area where pressure is more easily tolerated, and less pressure in the lower thoracic area where there is very little natural body padding and where even slight pressure is difficult to accept. (Refer to the preceding article by Campbell.)

Dr. Paul Campbell of the Portland Orthopedic Clinic was favorably impressed with Dr. Stuttle's presentation, and upon his return to Portland from the Academy meeting took immediate steps to make this orthosis available for trial.

A working model of the Stuttle orthosis was obtained from Plattner

Orthopedic Appliance Company of Peoria, Illinois, the facility which makes Dr. Stuttle's orthoses.

Since 1962 the Portland Orthopedic Clinic's Orthotic Facility has produced several hundred of these orthoses with great satisfaction to the patient, gratifying results to the orthopedic surgeon, and few problems to the orthotist.

## MATERIAL

The basic metal used is 20-24-T4 aluminum. The pelvic and thoracic bands are cut from sheet stock .080 thick. The two lateral uprights which extend from the end of the pelvic band upward to within half the width of the thoracic band are one-eighth inch by three-fourths inch 20-24-T4 aluminum bar. The two posterior uprights are attached to the lateral uprights three inches from where the lateral uprights are attached to the pelvic band.

---

<sup>1</sup> Chief Orthotist, Portland Orthopedic Clinic, Portland, Oregon.




Figure 1  
Measurements Required for Fabrication of the Stuttle orthosis.

### FABRICATION PROCEDURES

Procedures for fabrication of Stuttle orthosis are given below.

Measurements needed are:

1. Circumference around hips at the level of the greater trochanters (Fig. 1, left).
2. Circumference around waist just above crest of ilium (Fig. 1, center).

3. Length of thoracic band—this band should be only long enough to bear pressure on the flat surface of the back. It should not extend around the rib cage (Fig. 1, right).

The pelvic band (Fig. 2, left) as a rule will vary in length between twenty inches for a very large person to sixteen or seventeen inches for smaller individuals. The length of the




Figure 2  
Left, the pelvic band; center, the thoracic band; right, determination of the length of the orthosis.



pelvic band is easily determined. It extends around the hips at the greater trochanter level, but must not touch or exert any pressure on the greater trochanters. Straps, of course, do extend around and over the area of the greater trochanters without any problems, but the pelvic band must not exert any pressure on the greater trochanters.

The thoracic band (Fig. 2, center) is never longer than 11 in., the average being between 9 and 10 in. Iron wagon head rivets 3/16 in. dia. are used to fasten the joints, one rivet at each joint, in such a manner as to provide free motion.

Length of the orthosis is determined by measuring from the lower side of the pelvic band to the upper side of the thoracic band. The thoracic bands should always extend to the T11 level and should, if possible, extend to a point between the T9 and T10 levels (Fig. 2, right).

It is important to note that the lateral uprights are attached only to the pelvic band, allowing the top of the upright to swing free. Heat is needed to bend the posterior uprights away from the lateral uprights and the 45° bend to place the posterior uprights on an upward plane to be attached two and one-half inches from each end of the thoracic band (Fig. 3).

Holes for straps and "D" rings are placed at each end of the pelvic and thoracic bands and at the upper end of the lateral uprights.

Heavy leather straps are attached at the pelvic band, extending forward, just behind the point of attachment of the lateral upright.

Another heavy strap is attached at the top of the lateral uprights ex-




Figure 3

View of the Stuttle orthosis during fabrication to show configuration of the uprights.

tending back towards the end of the thoracic band where "D" rings are placed. The strap is threaded through the "D" rings and then brought around to the front where the top straps and the lower straps are attached to a heavy leather apron by four buckles placed in a manner that makes attachment easy. Average dimensions for the leather apron are 9 x 4½ in., and the buckles are placed more or less at each corner.

All that remains is covering, padding, and fabrication of the elastic apron. We all know it helps to have a nicely finished orthosis. It is some-




Figure 4  
The Completed Stuttle Orthosis

times difficult for patients to accept the idea of wearing a support, but if the orthosis is finished attractively, the patient is much more apt to accept it.

We cover the outside of the orthosis with white elk leather. The inside is lined with soft pearl cowhide or horsehide. Most orthotists have their own ideas on what they like to use for padding, and if they have a material that they like to use, it certainly could be used on this orthosis. We use one-half inch thick sponge rubber for padding the thoracic band, the pelvic band, and the lateral uprights. It is important to remember when covering the lateral uprights to leave a one-inch flap of leather extending forward for attachment of the elastic apron.

#### THE ELASTIC APRON

The elastic apron is made from Swiss weave elastic, eight inches wide with an eight inch open detachable zipper for a closure.

To determine the length of each side of the apron, the pelvic and waist measurements are needed.

The formula is:

From the circumference of hip, subtract six and divide by four to give one-half the length of the apron at the hip.

The same procedure using the waist measurement gives one-half the length of the top of the apron.

When the waist or hip circumference is uneven, *subtract* one inch to make it an even number.

For example:

The waist measurement of a	
patient is	29"
	and hips 37"
Subtract one inch from the	
37", making it	36"
(Only if an uneven number;	
if even number do not sub-	
tract)	Then subtract 6"

Then divide by 4" 30"

This gives you the lower



length of one side of your  
apron.  $7\frac{1}{2}"$

Now, do the same for the waist.

It being an uneven number,  
subtract one inch  $28"$

Subtract  $6"$

---

Divide by 4"  $22"$

$5\frac{1}{2}"$

So your apron would consist of  
two pieces of Swiss weave elastic:

$5\frac{1}{2}"$  at the top

$7\frac{1}{2}"$  at the bottom

An eight inch detachable zipper  
is sewn to one end of each piece  
and each piece is sewn to the flap on  
the lateral uprights. The complete  
orthosis is shown in Figure 4.

Once a general knowledge of this  
construction of the Stuttle orthosis  
has been acquired, much satisfac-  
tion can be gained by the orthotist.  
It is our hope that you will find this  
orthosis as gratifying to produce as  
we have at the Portland Orthopedic  
Clinic's Orthotic Department.

# Polyethelene Girdle for the Milwaukee Brace

By Jerome E. Skahan, C.O.,<sup>1</sup> John J. Skahan, C.O.,<sup>1</sup>  
and Walter E. Dase, Technician <sup>1</sup>

Since the inception of the "Milwaukee Brace" orthotists have been searching for a material to substitute for the leather originally recommended for fabrication of the girdle.

Leather can be molded readily, it can breathe, and is quite compatible with the outer body tissues. However, leather has its drawbacks. Good grades are much in demand, and not always readily available. An otherwise good hide may have small imperfections which make it unusable. Leather must be mulled and molded while wet. The drying time is long, and molded leather must be lined with horsehide or some equivalent. The procedures are time consuming, and therefore expensive.

Some plastics with promising characteristics were presented a few years ago for replacement of leather. Unfortunately, these materials had

drawbacks that were not noticed immediately. Experiments with these plastics in our facility showed them to be acceptable but not ideal.

A material that we have found to be most satisfactory is low-density Polyethelene. Polyethelene has a smooth finish, and it has presented no problems of chafing or reddening the skin owing to abrasion, especially in the waist area. It does not need to be lined. It has the "snap" and flexibility similar to leather so that it can be bent, yet does not lose its ability to control. More than 30 braces using Polyethelene have been fitted during the past year.

Polyethelene is commercially available in 4' x 8' sheets, 3/16" thick, at a cost of less than \$20.00 per sheet. When properly laid out at least ten pieces 12" x 36" can be cut from one sheet and there will still remain material for making pads, etc. (Fig. 1).

<sup>1</sup>Central Orthopedic, 2031 Auburn, Ave., Cincinnati, Ohio 45219.



$36'' \times 12''$			$36'' \times 12''$			$24'' \times 12''$	
$36'' \times 12''$	$36'' \times 12''$	$36'' \times 12''$	$36'' \times 12''$	$36'' \times 12''$	$36'' \times 12''$	$36'' \times 12''$	$36'' \times 12''$

Figure 1

Layout for 4 ft. by 8 ft. sheet of Polyethylene sheet to obtain maximum number of girdles for the Milwaukee Brace.

### TOOLS AND MATERIALS NEEDED

1. Forced air oven capable of 230° F.
2. Cloth covered wood frame.
3. Commercial Talc (Merck).
4. Heavy elastic strap.
5. China marker.
6. 2-3 pr. insulated gloves.
7. Hammer and nails.
8. Knife.
9. Heat gun.
10. Felt cone or muslin wheel.

The cast is smoothed and vertical-fabrication lines are marked anteriorly and posteriorly. Trim lines can be marked on the cast also at this time (Fig. 2). A piece of the plastic is cut to the desired height and length according to the measurements of the cast. A vertical centerline is marked on the plastic with

a china marker or similar marking pencil. The plastic is placed in an oven pre-heated to 230° F, for eight minutes (Fig. 3). Use a well-powdered heavy cloth stretched over a wood frame for a non-adhering surface. After the plastic has been heated properly it is carefully removed from the oven and transferred to the cast, making certain that the center line on the plastic is aligned with the anterior vertical fabrication line on the cast. An elastic strap, 2" wide and 36" long, is stretched around the waist while the rest of the plastic is held manually next to the cast. Nails are used along the posterior edges to hold the plastic to the cast (Fig. 4).

*N.B. Avoid overheating.* If the material is allowed to heat for a longer period of time or at a higher heat than indicated, it will distort e.g., reduce in length and width but




Figure 2  
Trim lines marked on cast for the Milwaukee Brace.

become thicker, softer than desirable, tacky, and the surface will become uneven or bubbly. While the material is still quite warm, the excesses posteriorly can easily be trimmed with a knife.

Three people are needed for a very short period for the actual molding procedure. One person is needed on each side of the cast in order to position the hot plastic while a third person stretches the elastic strap into the waist area. This same person can then hand the elongated ends of the strap to the people on each side while he nails the plastic to the cast. The elastic strap is then removed and the excess plastic material is cut away.

After it has cooled and "set," the

plastic may be cut along the trim lines marked on the cast (Fig. 5) with a heat gun and sharp knife. Because Polyethylene is translucent when heated the lines are visible at this time. The edges are smoothed initially on a felt cone and, finally, on a muslin wheel at 3450 RPM. The material will not "burn" or melt at this speed.

The attachment of anterior and posterior uprights to the Polyethylene girdle is done with jiffy, or speed rivets (Fig. 6). The use of copper rivets for attachment is not recommended, because the plastic can be cut slightly by the head of a copper rivet, thus inducing splitting of the plastic.

It has not been necessary to use






Figure 3  
Polyethylene sheet being put into oven on a cloth-covered wood frame.




Figure 4  
Use of elastic strap for molding the Polyethylene sheet, left, and anterior view showing the use of nails to hold the sheet in proper position over the mold.




Figure 5  
The Polyethylene girdle after being trimmed on the cast.




Figure 6  
The Anterior and Posterior Uprights in Place.

a pelvic hinge on any of the Milwaukee braces made with Polyethelene. There have been no fractures or failures due to repeated flexure.

Polyethelene is washable and is impervious to most cleaning solutions. Heating with a heat gun permits slight adjustments. Polyethelene of this same thickness can also be used for making the correctional pads for scoliosis and kyphosis. The pads are lined with Plastazote and a Dacron strap is used for attachment.

Plastazote is washable also.

Low-density Polyethelene has answered many of our problems concerned with the girdle of the Milwaukee brace. It is inexpensive, relatively easy to form, easy to keep clean, relatively easy to adjust, and is flexible but firm.

Polyethelene is not considered to be the last word in materials for the girdle of the Milwaukee brace but at this time appears to be the best for this purpose.




# **Use of a Cast-Brace Application Device, Plaster Brim, and Ankle Support for Fracture Bracing**

Max Lerman, C.O.<sup>1</sup>


A device has been designed to permit one man to apply cast braces for lower-limb fractures (1, 2) in a faster time and with more accuracy than is the case with present procedures.

---

<sup>1</sup> Lerman and Son, 8710 Wilshire Boulevard, Beverly Hills, Calif. 90211.





The present joint-alignment fixture is used as the basic knee support. When the joint alignment fixture is reversed it can hold the knee as in a vise. It can be raised or lowered to accommodate some knee flexion. A tilt adjustment feature makes it possible to accommodate for knee or tibial rotation. An ordinary crutch handgrip serves as padding under the popliteal area.




Stockinette is applied over the usual interface material, and the foot is suspended by tying the stockinette to the suspension bar. When the heel is centered with the device the knee hinges will always be horizontal. When valgus or varus knee stress is required, the stockinette may be attached to an outrigger, which can be reversed so as to serve on each side.



A plaster brim (see below), a Truform No. 782 cloth ankle support lined with  $\frac{1}{4}$  foam rubber, and relief padding are applied.



The proximal and distal casts are applied.



Knee hinges are applied outside the device. A clamp with screws above and below holds the hinges firmly in place. The hinges are set slightly posterior to midline. The use of hose clamps is recommended. The thickness of the knee jig uprights provides the proper clearance between the knee hinges and the knee.



The hinges are secured with plaster in the usual manner.

I find that the cast-brace alignment device helps us to make the brace better fitting and more comfortable. It also saves time by making a difficult job easier.

The development of the cast-brace application device would not have been possible without the cooperation of Dr. J. P. Harvey of the Los Angeles County Medical Center, University of Southern California and Dr. Vert Mooney of Rancho Los Amigos Hospital, under whose auspices and guidance cast braces are applied.


U.S. Manufacturing Company, Glendale, California will soon distribute this new cast-brace application device.





### **THE PLASTER BRIM**

The plaster brim is made out of one roll of 4 inch or 6 inch plaster over stockinette using a plastic brim for the quadrilateral shape. The opening is on the lateral side. Use of plaster has reduced skin problems because of its breathing properties, and made trimming and relieving easier.



### ANKLE SUPPORT

When the ankle support is used problems of ankle-joint alignment are avoided because slight motion in pronation and supination, plantar and dorsiflexion are allowed. The foot is in normal contact with the floor while weight-bearing and an ordinary slipper can be worn. Heel pressure sores have been eliminated. The brace weighs about three pounds less. The use of the ankle brace is limited to femoral, knee and tibial plateau fractures. It is not indicated in the case of fractures located more distally.\*

### LITERATURE CITED

1. Committee on Prosthetics Research and Development, Cast-bracing of fractures National Academy of Sciences, Washington, D.C. 1971.
2. Mooney, V., V. L. Nickel, J. P. Harvey, Jr., and R. Snelson, Cast-brace treatment for fractures of the distal part of the femur. *J. Bone Joint Surg.*, 52-A:8:1563-1578, December 1970.

\* It should be used in selected cases only, and is by no means a substitute for the rigid plaster. I lace the anklet proximal to distal so that the patient can adjust the laces when he is reclining.

# Use of Photography in Cast Preparation

James H. Tyo<sup>1</sup>  
George Hall, Jr., C.P.<sup>2</sup>

Cast modification may well be considered the single most important phase of prosthetic service because it is the key to the provision of a comfortable, functional socket. This is, obviously, the basis for all prosthetic treatment.

The prosthetist draws on a wealth of knowledge in preparing the plaster model, not the least of which is his memory of each individual amputation. It was felt that a method of recording the visual appearance of the stump at the time of casting might be quite helpful to the prosthetist.

The Polaroid "Big Shot" camera (Fig. 1), type 108 film, and battery-less magicubes were selected because of low initial investment, con-

venient operation, and the immediate results provided by the Polaroid process. The time added to the normal casting and measurement procedure is thus less than four minutes.

The basic procedure is quite simple. When possible, we position the patient against a dark background, normally a piece of brown paper taped to the cast room wall, and place an 18 in. ruler parallel to and directly beside the amputated side. On the stump, we mark the reference point for measurement of length (i.e., on the BK the patella tendon). The photo is taken from the front (Fig. 2).

There are obvious advantages when the photographs are used, the visual reminder of shape and condition being the most important (Fig. 3). There is, however, a secondary advantage. During the testing of this procedure, we noted an immediate

<sup>1</sup> Prosthetic Services, Frees and Tyo, Inc., 1124 East Fayette St., Syracuse, New York 13210.

<sup>2</sup> Prosthetic and Orthotic Services of Buffalo, C. B. Winn, Inc., 699 Hertel Avenue, Buffalo, N.Y. 14207.





Figure 1  
The 'Big Shot' camera.




Figure 2  
The frame shown attached to the camera has been designed to eliminate the problem of focusing at an uncomfortable level.



Figure 3

The photograph is attached to the patient's file as a 'visual' record of the stump. Its primary use is during the cast modification procedure.

improvement in patient response to treatment.

Our feeling at this time is that the patient responds to the quickly available photograph as a demonstration of our interest in offering the best prosthetic care possible. This, of course, is a highly subjective analysis that can only be proved by the passage of time.

This addition to the normal cast and measurement procedure should be viewed as just that, an addition.

It is not a replacement, or substitute for thorough prosthetic examination of the stump, and not a "short-cut" method of recording prosthetic information. The procedure has been developed as an accurate, simple, and inexpensive means of providing a picture of the amputee's stump while preparing the plaster model. The prosthetist can "see" his patient, thus adding to the written information and his memory of shape, structure and texture.

# Population and Access to Prosthetists<sup>1</sup>

G. E. (Ned) Snarples, Ph.D.<sup>2</sup>

Concern with delivery systems for health care has become widespread within government and within many of the professional groups involved in the provision of care. As part of this concern, much attention has been focused upon locating health facilities so as to enable maximum access for populations in need of services. It is recognized that location factors do not produce or explain utilization, nor necessarily effect quality of health services. Still, the location of facilities is a basic determinant of *potential* access to

services, whether or not utilization occurs. Without a facility to use there is no utilization. With a facility accessible, there may be some utilization. Judgments of the adequacy of services begin, but do not end, with the question of availability or potential success within geographic territories.

This paper describes the geographic dispersion of certified prosthetists and certified facilities in the United States. It also presents an imperfect but useful analysis of the potential access which certain general populations have to prosthetists, based upon ratios of facilities to population. In general, the accessibility of prosthetic services seems extremely uneven.

Caution in interpreting these findings must be exercised. Again, results do not necessarily relate to quality of services and no such implication should be drawn. Also, there are

---

<sup>1</sup> Based on research supported by the Maternal and Child Health Service under grants PC-1003 and MC-R-260044-05-0 titled: Child Amputees: Disability Outcomes and Antecedences.

Major contributions to this analysis were made by Miss Rochelle Coonley and Miss Roseanne Allevato. Their help is gratefully acknowledged.

<sup>2</sup> Assistant Professor of Maternal and Child Health and Director, Child Amputees Research Projects, School of Public Health, The University of Michigan.



many additional factors which influence actual accessibility and real utilization of prosthetic services. Among these are the additional number of competent but uncertified prosthetists, the restriction of practice by some prosthetists to institutionalized or other special populations (e.g., VA), transportation, economic, and many other social factors influencing decisions to use health services.

While this analysis cannot yield conclusions about the quality of services nor the efficacy of utilization in relation to population needs, it is useful in several ways. It identifies a number of geographic areas and cities where there is reason to believe that prosthetic services are insufficient by being absent, or nearly so. Questions about service quality in such areas are unnecessary if there are no services. The analysis also identifies areas and cities where further investigation of actual accessibility would be likely to be most fruitful. This preliminary work provides a basis for selecting areas or cities for comparison purposes in further work seeking to define adequate prosthetic services for populations. It offers a starting point to the profession for self-examination and eventually may enable the profession to set goals, plan education programs, and influence the availability of services.

Our interest in the access problem results from intensive interviews with upper- and lower-extremity amputees throughout the country during the last five years. The more than 500 amputees with whom we talked gave clear evidence that the accessibility of adequate prosthetic

services has considerable influence upon continuity of prosthesis wearing, decision to permanently discontinue use, and general levels of physical and social performance. Further, we heard claims that access to prosthetic services varies considerably by geographic areas.<sup>3</sup>

To better comprehend the availability of such services, we plotted locations of all certified limb facilities on a map of the United States. This exercise clearly demonstrated that tight clusters of facilities occur in and around certain big cities and that many large areas of the country have no certified prosthetists. Such a pattern can be expected because of the variation in population density.

The question then becomes: Does the distribution of prosthetic services match the distribution of population such that potential access is equal throughout the country? A further question, not investigated here, is: Are the human problems of obtaining prosthetic services within reasonable limits throughout the country?

The second question requires original field data which can only be obtained at considerable expense, requiring special funding. The first question can be attacked using al-

---

<sup>3</sup> Experiences of the amputees we interviewed were not known to be representative of all amputees, or of those in any one area. Our respondents were the widely dispersed patients of a single clinic caseload. Their reports could not yield a systematic picture of the availability of prosthetic services within any one geographic area, or the country as a whole. But they were highly suggestive. Patients were dispersed from Alaska to Puerto Rico, from Los Angeles to New Hampshire, from North Dakota to Texas. In all, 345 people were extensively interviewed, less extensive data was collected for over 500.

ready existing information. The issue is whether the geographic distribution of facilities is consistent with the distribution of people needing prosthetic services. An estimate of this relationship can be made by comparing population figures to the number of prosthetists in the same area.<sup>4</sup> There is certainly no ratio of population to prosthetists which is known to be satisfactory; that is, representative of optimum service accessibility. It would be hazardous to estimate the size of population which could be satisfactorily served per average prosthetist.

However, we can still detect regions with high and with low ratios of population to certified prosthetists and interpret the ratios using added knowledge of areas and services.<sup>5</sup> Variation in a series of such ratios for different areas would suggest, at the least, that in some areas prosthetic services are potentially less available than in others. Extreme variation would suggest that inadequate availability exists in some areas, if need is constant, and/or oversupply might exist in others, the latter being much less likely according to our respondents. Using existing data, and avoiding the problem

of establishing an expected patient volume per prosthetist, our method of comparing ratios is as close as we can practically and inexpensively come to estimating whether there are satisfactory numbers of prosthetists for populations.

Following this logic we sought to learn whether access to prosthetists is likely to be a problem in the 25 most populous metropolitan areas of the country. United States Census population figures for 1970 were obtained for each metropolitan area.<sup>6</sup> Counts were made of all certified prosthetists and facilities in each metropolitan area using the *1971 Registry for Certification* published by the Board for Certification. Ratios of total population to certified prosthetists and ratios of total population to certified facilities were computed for each area and are shown in Table I.

Inspection of the ratios reveals a striking degree of variation even among our most populous areas. Within our 25 largest cities, the highest ratio was almost eight times the lowest for the number of facilities. The number of certified prosthetists varied by a factor of about three and one-half. Thus, Newark has eight times as many people per certified facility as Cleveland. Milwaukee has about three and one-half times as many people per certified prosthetist as Cleveland.

<sup>4</sup> Provided we accept certain assumptions specified later.

<sup>5</sup> For example, certain facilities utilize personnel in ways that effectively serve more clients, especially using uncertified prosthetists. The numbers of such people may vary in different regions. Cities vary in the degree to which their services are depended upon by people living outside the metropolitan area limits defined by the census. Reputations about quality and speed of services are certainly relevant qualifications to any ratios. I am indebted to A. Bennett Wilson, Jr., Executive Director, Committee on Prosthetics Research and Development, for examples illustrating this important point. The need for a thorough study by AOPA should be obvious.

<sup>6</sup> Population figures were obtained from the U.S. Bureau of Census, *Number of Inhabitants* (1970), by state. The AOPA 1971 listing of individual prosthetists was used. Hand counting was done carefully but may be subject to slight error. Census data for each state was used to compile lists of political divisions included in each SMSA. Judgment was required when a listed address was near an SMSA which was incompletely defined.



TABLE 1

# POPULATION RELATIVE TO CERTIFIED PROSTHETISTS AND FACILITIES IN THE 25 LARGEST U.S. METROPOLITAN AREAS

Area	STANDARD METROPOLITAN STATISTICAL AREAS	CERTIFIED		RATIOS	
		1970 Total Population	Facilities	Pros- thetists	Population in Thousands Facilities Prosthetists
New York	11,571,899	12	49	964:1	236:1
Los Angeles-Long Beach	7,032,075	17	25	414:1	281:1
Chicago	6,978,947	10	23	698:1	303:1
Philadelphia	4,817,914	8	13	602:1	371:1
Detroit	4,199,931	4	12	1050:1	350:1
San Francisco-Oakland	3,109,519	8	13	389:1	239:1
Washington, D.C.	2,861,123	3	10	954:1	286:1
Boston	2,753,700	6	11	459:1	250:1
Pittsburgh	2,401,245	2	5	1201:1	480:1
St. Louis	2,363,017	3	9	788:1	263:1
Baltimore	2,070,670	3	8	690:1	259:1
Cleveland	2,064,194	1	3	2064:1	688:1
Houston	1,985,031	2	5	993:1	397:1
Newark	1,856,556	7	7	265:1	265:1
Minneapolis	1,813,647	5	9	363:1	202:1
Dallas	1,555,950	2	5	778:1	311:1
Seattle-Everett	1,421,869	3	6	474:1	237:1
Anaheim-Santa Ana- Garden Grove	1,420,386	0	4	—	355:1
Milwaukee	1,403,688	3	7	468:1	201:1
Atlanta	1,390,164	2	5	695:1	278:1
Cincinnati	1,384,851	2	3	693:1	462:1
Paterson-Clifton-Passaic	1,358,794	1	3	1359:1	453:1
San Diego	1,357,854	2	3	679:1	453:1
Buffalo	1,349,211	1	2	1349:1	675:1
Kansas City	1,253,916	2	3	627:1	418:1

SOURCES: Population figures were obtained from the U.S. Bureau of Census, *Number of Inhabitants (1970)*, by state, 1971 *Registry of Certified Facilities and Individuals in Orthotics and Prosthetics*, American Board for Certification in Orthotics and Prosthetics, Inc., Washington, D.C., 1971.

The size of the population per prosthetist is large. Table 1 ratios are in thousands of people for each certified prosthetist. Thus, Cleveland has 688,000 people for each certified prosthetist. While the presence of apprentices and uncertified prosthetists, some working in certified facilities, would reduce the ratios, it remains a rough index to service access.

Inspection of Table 1 reveals some differences concomitant with population size among the metropolitan areas. The most populous areas (the first eight listed, New York through Boston) have predominantly middle ratios (five) and

low ratios (three), presumably having average or better access to prosthetists, among the 25 areas compared. The eight metropolitan areas which fall in a middle category according to population size (Pittsburgh through Dallas on the list in Table 1) are mixed in ratios, low (four), middle (two) and high (two), generally having average or better access to prosthetists. The nine metropolitan areas with the smallest populations among these 25 (Anaheim through Kansas City) have predominantly high ratios; five of the nine are high ratio cities, two each are low and middle ratio. These ratios may indicate less access



to prosthetists.

Unlike any other consecutively listed metropolitan areas in the list, the last five areas all fall at one extreme and all have high population ratios. While this pattern is weak, it may indicate a trend towards higher ratios and less access as city size decreases, perhaps after a threshold level of population is reached.

It is reasonable to expect that even greater variation would be revealed if ratios were computed for smaller cities. This would in part be due to the small numbers of prosthetists but also to the absence of services in many populous areas.

This does not suggest that any specific ratio reflects adequate access. It is most probable that no city has an adequate supply of prosthetists in relation to the numbers of people with deficiencies and the services they need. This is what our interviews suggest.

It must be noted that this method of estimating the relative availability of prosthetic services rests upon the following assumptions. The total number of people in an area is used as an index to the number in that population with prosthetic service needs. Comparisons between cities are based on the assumption that, among our largest metropolitan areas, the proportion of each population which is afflicted is probably about the same. It is also assumed that because there are concentrations of facilities and prosthetists in these areas, the differences between them in the number of people who are served per prosthetist will average out. If it does not average out, then there is chance of misin-

terpretation of the ratio variation, provided that those cities with high populations per prosthetist are the same cities in which prosthetists are able to serve more people.<sup>7</sup> It seems unlikely that there are sufficient differences in the numbers of patients served per average prosthetist and that these differences occur in a pattern to offset the effect of population ratios.

If these assumptions can be accepted as reasonable, then a nationwide look at the variation between ratios is warranted. Table 2 groups cities roughly into the one-third having the lowest ratios, the one-third having the highest ratios, and a middle one-third. Those with the lowest ratios, having the smallest population for each prosthetist, are likely to be cities where prosthetists are more accessible, based on sheer numbers alone. Metropolitan areas with the highest ratios probably have less accessible prosthetic services. It is interesting that no region of the country can be identified, from Table 1, as having either mostly cities with relatively easy access to prosthetists or as having mostly cities with rela-

<sup>7</sup> If the number of people who are served by low volume prosthetists differ from the number served by high volume prosthetists by a factor of three and one-half or more, then we must modify conclusions based on the finding of variations by a factor of three and one-half in the ratios of prosthetists to population. As an example, if a high population ratio city such as Cleveland was known to have high volume prosthetists serving three and one-half times as many patients as a low population ratio city such as Milwaukee, where it was known that low volume prosthetists were serving less than one-third as many patients as were served by the Cleveland prosthetists, then the three and one-half fold difference in prosthetist to population ratios would be offset or equalized and a conclusion of differing availability would be inaccurate.

**TABLE 2**  
**GROUPED DATA—CERTIFIED PROSTHETISTS**  
**RATIOS TO POPULATION**

**Low Ratios—Greatest Access—201:1 to 265:1 (thousands)**

Milwaukee	Seattle	Baltimore
Minneapolis	San Francisco	St. Louis
New York	Boston	Newark

**Middle Ratios—Average Access—278:1 to 371:1 (thousands)**

Atlanta	Chicago	Anaheim
Los Angeles	Dallas	Philadelphia
Washington, D.C.	Detroit	

**Highest Ratios—Least Access—397:1 to 688:1 (thousands)**

Houston	Paterson	Buffalo
Kansas City	Cincinnati	Cleveland
San Diego	Pittsburgh	

SOURCES: Population figures were obtained from the U.S. Bureau of Census, **Number of Inhabitants (1970)**, by state. **1971 Registry of Certified Facilities and Individuals in Orthotics and Prosthetics**, American Board for Certification in Orthotics and Prosthetics, Inc., Washington, D.C., 1971.

tively difficult access. A mixture of both high and low ratio metropolitan areas appears in each general zone.

A listing of cities and ratios for each AOPA region, Table 3, reveals the same mixed results in general. The exception is that, in Region V, two of the three large metropolitan

areas have high population ratios, Cleveland (688:1) and Cincinnati (462:1), and the third has only a middling ratio, Detroit (350:1). An opposite pattern occurs in Region VI where two low ratio cities, Milwaukee (201:1) and St. Louis (263:1) occur with a middle range

**TABLE 3**  
**POPULATION PER PROSTHETIST RATIOS**  
**BY REGIONS OF AOPA**

I Boston	250:1	L	VI Milwaukee	201:1	L
			Chicago	303:1	
II New York	236:1	L	St. Louis	263:1	L
Buffalo	675:1	H			
Newark	265:1	L	VII Minneapolis	202:1	L
Paterson	453:1	H	Kansas City	418:1	H
III Philadelphia	371:1		VIII Houston	397:1	H
Washington, D.C.	286:1		Dallas	311:1	
Pittsburgh	480:1	H			
Baltimore	259:1	L	IX Los Angeles	281:1	
IV Atlanta	278:1		Anaheim	355:1	
			San Diego	453:1	H
V Detroit	350:1		X San Francisco	239:1	L
Cleveland	688:1	H			
Cincinnati	462:1	H	XI Seattle	237:1	L

KEY: H = High ratio cities, having larger populations per prosthetist (397:1 to 688:1)—thousands

L = Low ratio cities, having smaller populations per prosthetist (201:1 to 265:1)—thousands

SOURCES: Population figures were obtained from the U.S. Bureau of Census, **Number of Inhabitants (1970)**, by state. **1971 Registry of Certified Facilities and Individuals in Orthotics and Prosthetics**, American Board for Certification in Orthotics and Prosthetics, Inc., Washington, D.C., 1971.



city, Chicago (303:1). In the remaining regions where there are two or more metropolitan areas, high population ratio cities appear with offsetting low ratio areas. Middle range areas are well scattered among all regions.

Combining the ranking of ratios for prosthetists and for facilities we can list metropolitan areas in which there may be scarcity of both shops and certified prosthetists. Note that this conclusion is tentative, at best, and subject to modification by detailed knowledge of the quantity and character of services. The potentially difficult access areas are Cleveland, which ranks at the extreme on both prosthetist and facility ratios, Buffalo, Pittsburgh, Houston, and Paterson. Paterson is the only one of these geographically near another area with ratios likely to be offsetting. The cities which are likely to have the readiest access, with low population ratios for both individual prosthetists and facilities, are Seattle, Milwaukee, Boston, San Francisco, Minneapolis, and Newark. Again, Newark may be used by Paterson area clients, reducing its real access.

It appears from this analysis that, even among populous metropolitan areas, we find sufficient variation in the numbers of certified prosthetists per population to expect important differential effects upon the availability of services and, therefore, upon the well-being of patients. It seems likely that the variation in access is even more extreme among smaller metropolitan populations. It is important not only to document this distribution more fully but to learn in more direct and concrete terms

about the problems of access to services.

Returning to our basic question about the reality and effects of thinly spread prosthetic facilities away from the city clusters, we must re-emphasize the need for research on the human problems of obtaining services. From the analysis of metropolitan areas it seems clear that the problems of obtaining access to prosthetists are not limited to people living in ranch and farm areas. While research on personal problems of access is more complicated and must await the mounting of considerable research resources, it is useful to outline the dimensions of the problem for areas where the number of prosthetists is low.

Half of our states have only three or fewer facilities certified in prosthetics; twenty-four states have only two or less, fourteen have only one, and seven states have no facility certified in prosthetics. If both prosthetic and orthotic facilities are counted, the picture is only slightly better. Six states have no facilities, of either kind, thirteen have one or less, and seventeen states have two or less. Lists of these states are shown in Tables 4 and 5.

Generally, the states with few facilities are relatively less populous, either relating to low density or small size. However, many are likely to have a higher proportion of amputees in the general population, based on the contributions of certain etiologies to the total afflicted. The long-standing leading causes of amputation have been farm machine accidents and gunshot (primarily hunting) accidents. While these causes are dropping in rank among



**TABLE 4**  
**CERTIFIED PROSTHETIC FACILITIES**  
**IN STATES WITH FEW FACILITIES**

(Excludes Shops for Orthotics Only)

No Facilities	One Facility	Two Facilities	Three Facilities
Alaska Maine New Hampshire North Dakota Rhode Island South Dakota Wyoming	Delaware Hawaii Idaho Kansas Nebraska Nevada Vermont	Alabama Arizona Arkansas District of Columbia Mississippi Montana New Mexico Oregon Utah West Virginia	Kentucky Oklahoma
Total = 7 States	Total = 7 States	Total = 10 States	Total = 2 States

NOTE: The number of certified facilities per state indicates whether or not there is any possibility of in-state quality service for a population of whatever size. It is especially important because of the changing sources of payment and referral from or through state agencies.

**TABLE 5**  
**CERTIFIED PROSTHETIC AND/OR ORTHOTIC FACILITIES**  
**IN STATES WITH FEW FACILITIES**

No Facilities	One Facility	Two Facilities
Maine New Hampshire North Dakota South Dakota Wyoming Rhode Island	Nevada Idaho Nebraska Vermont Delaware Alaska Hawaii	Montana Arizona New Mexico West Virginia
Total = 6 States	Total = 7 States	Total = 4 States

reasons for new amputations, they probably are still the causes which contributed the most cases among all existing (prevalent) cases. These causes occur mostly in rural areas and therefore we might expect a higher proportion of amputees in rural, less populous, states. In addition, these same states tend to have older populations and therefore may have more geriatric amputees, relative to total population, than states having more prosthetic facilities.

The lack of facilities in some states is to some degree compensated

for by the existence of other facilities nearby. For example, the low number of facilities in Rhode Island, Delaware, and the District of Columbia must be viewed in the perspective that concentration of facilities occur in Boston and Hartford, Philadelphia, and Baltimore. The question of whether the human problems of time, cost and inconvenience deter access for people within these densely populated areas must await specific study.

Many states without suitable prosthetic facilities are not adjoined by

other states having a concentration. These become zones of sparse prosthetic resources. In New England, three states of the six have no certified prosthetic facility and a fourth (Vermont) has one, located in the extreme northwest corner of the state. In effect, all of Maine and New Hampshire and two-thirds of Vermont must be served from Boston. Most dramatic, all of both Dakotas, half of Minnesota, half of Montana, most of Nebraska, all of Wyoming, and much of Kansas, all geographically contiguous, are without a certified prosthetic or orthotic facility. It seems doubtful that the scattered shops next to this vast area can offer adequate access to services. This leaves a zone in which there are no certified facilities that is about 1,000 miles wide east to west at the top and 1,000 miles long from north to south. The concentrations of cer-

tified facilities nearest to this giant vacant V shape in our northern mid-west occur at Chicago, Denver, Oklahoma City and Minneapolis.

Population figures and demographic descriptions such as those used here are not direct measures of prosthetic service accessibility. However, they provide means to estimate the existence, extent, and location of problems in access. It is our conclusion that access to prosthetic services is clearly unequal and deserves national attention. Adding what is known about the effectiveness of shops and individual prosthetists, professional judgments can be made defining the problems and establishing their priorities. Hopefully, resources will be allocated not only for the research needed to describe the human problems of access to prosthetic services but also to alleviate those problems.


**A FEW**

## ***Fillauer Favorites***

**FILLAUER...The Specialists in Specialties**

**FROM**

**TRILATERAL SOCKET  
PERTHES ORTHOSIS COMPONENTS**




**All Immediately Available—Quality At Low Cost**

**TO**

**PLASTIC LAMINATE ISCHIAL  
WEIGHT BEARING CUFFS**

**Sizes 2-4½"—Order #3302**

**Specify Size AP, Right or Left**



**DISPOSABLE GLOVES TO  
PROTECT YOUR HAND**

**ALL FROM**

***Fillauer***

**AND EVEN**



**Write for New "Green"  
Catalog**

**Established 1914**

**Box 1678**

**Chattanooga, Tennessee**





**KNIT-RITE, INC.**

1121 GRAND AVENUE • KANSAS CITY, MISSOURI 64106

PHONE: 816-221-0206

## PLASTAZOTE SHOE



Removable inner sole may be easily heated and molded to the pressure, profile of the patient's foot, further improving pressure distribution for various foot deformities.

### PLASTAZOTE SHOE

This shoe was engineered in conjunction with leading physicians to fill a long standing need to relieve pressure in the diabetic, arthritic and other patients requiring special pressure relief. Immediate ambulation may be obtained in cases where Edema prohibits fitting of a standard shoe.

Available in six common sizes enabling a wide variety of feet to be fitted. The unique Plastazote material can be further molded or modified using a standard heat gun. When heat is applied in the required area the shoe can be reshaped by hand to provide relief for bunions, etc.


1. Excellent pressure distribution
2. Impervious to perspiration, common acids and alkalis
3. Perforated for ventilation
4. Washable
5. Immediate fitting
6. Comfortable
7. Economical
8. Lightweight
9. Available in sizes 6, 7, 8, 9, 10 and 11, Right or Left



# ORTHOLEN




## Posterior Leaf Spring Below Knee Orthosis



Designed to control most drop-foot conditions, the ORTHOLEN-Posterior-Leaf-Spring-B.K.-Orthosis features:

- Controlled dorsiflexion assist.
- Cosmesis: ORTHOLEN is a flesh-colored thermoplastic, fitting intimately against the back of the leg and a sock or a hose may be worn under as well as over the orthosis.
- Support and control of foot-architecture.
- Interchangeability of shoes.
- Light weight—5 sizes from 4 to 9 ounces!
- Ease of adjustability. ORTHOLEN can be readily thermo- or cold formed. Extremely solid.
- Hygienic.



The ORTHOLEN-Posterior Leaf Spring B.K. Orthosis is available in 5 sizes, left and right in a *pre-molded* form. It is recommended however that it be individually molded to a cast of the patient's leg.

- 1 — Premolded Form  
2 — Grinded and Fitted

*Distributed by:*

FILLAUER SURGICAL SUPPLIES INC.  
CHATTANOOGA, TENNESSEE 37401


R.T.C. INC., ORTHOTIC-PROSTHETIC COMPONENTS  
612 MONROE AVE., KENILWORTH, N.J. 07033

WILH. JUL. TEUFEL, NECKARSTRASSE 189-191,  
STUTTGART, W.GERMANY

MANUFACTURERS OF BANDAGES, ORTHOPEDIC AND PROSTHETIC ARTICLES



now...an external  
breast prosthesis  
so real it looks  
and feels like  
part of her.



There's no noticeable difference with **COMPANION**. An artificial breast made of soft, mobile silicone gel encased in a smooth silicone rubber skin — designed to match the normal human breast in contour, weight, fullness and feel.

Worn right next to the skin in an ordinary bra...the silicone rubber skin will actually take on the body temperature... this prosthesis feels like it's *really* part of the body.

For quality and competence that compliment your professional skills, prescribe the **COMPANION** breast prosthesis by **Truform** with confidence.

**TRUFORM**

*anatomical supports*

**CINCINNATI**

3960 Rosslyn Drive • Cincinnati, Ohio 45209 • (513/271-4594)

**SAN FRANCISCO**

589 Howard Street • San Francisco, California 95104 • (415/781-1135)

**NEW YORK**

4 West 16th St. at 5th Ave. • New York, N. Y. 10011 • (212/686-6460)

**IN CANADA**, Contact Airway Surgical Appliance Industries Ltd.,  
Ottawa, Ontario. K1R5T8



# Uptight orthopedic appliances

If your Truform orthopedic appliances are uptight, you won't be. These items are available in stock size ranges sufficiently wide that no difficulty will be encountered in finding an item that will fit perfectly.



Model No. 962

Chair-back brace, with full corset front, has no real equal for low-back disorders. The brace's metal framework of four-uprights and two horizontal bands has tremendous fixation potential.



Model No. 782

Cock-Up Wrist Splint is capable of delivering support to immobilize the wrist area. It is made of soft yet durable leather and comes in six sizes, either left or right. No. 771 is exactly the same, but the material is heavy canvas.




Model No. 9971-B

Hinged Knee Cap is now offered with Velcro fasteners. The cap's unique feature is that the application and adjustment are greatly simplified by the wide Velcro straps that close the appliance. The circular cutout over the knee eliminates irritating pressure in the back, and the cap comes with cartilage pads which may be removed if desired.



Model No. 975

Cervical Collar is called a universal collar because its two parts, joined together by means of Velcro closures, can be individually set in so many different positions that a collar of almost any shape or dimension can be made on the spot.



Model No. 1501

Economy Model Head Halter Home Traction Kit comes equipped with 1502 head halter and notched spreader bar, over-door traction hanger 989 complete with pulley, water weight bag 9960-S and a section of rope. Kit also available with 20 lb. water bag, if specified.

**Write: Truform Anatomical Supports, sold only by Ethical Appliance Dealers.**



TRUFORM RED BOOK

A complete catalog illustrating the entire Truform line.


## TRUFORM

anatomical supports

CINCINNATI  
5900 Resolyn Drive • Cincinnati, Ohio 45209 • (513/271-4394)  
SAN FRANCISCO  
589 Howard Street • San Francisco, California 94104 • (415/781-1135)  
NEW YORK  
4 West 18th St. 31 5th Ave. • New York, N. Y. 10011 • (212/686-6600)  
IN CANADA, Contact Airway Surgical Appliance Industries Ltd.,  
Ottawa, Ontario, K1R0T8

# Becker ADJUSTABLE HYPEREXTENSION BRACE

- Quick release, snap-out attachment
- Adjustable, self aligning posterior pad
- Rotating adjustment for sternal and pubic pads
- Vertical and horizontal sliding adjustments
- Bi-lateral worm gear traction bands
- Plastic water resistant pad covers
- Constructed of 24 ST aluminum



MODEL  
L-25

1776 South Woodward • Birmingham, Michigan

MANUFACTURERS OF PRECISION-MADE BRACE PARTS

BECKER ORTHOPEDIC APPLIANCE COMPANY

24 Hour Service

**Realastic®** The largest, most complete line of  
cosmetic restorations available anywhere in the world.

- COSMETIC GLOVES
  - PASSIVE & PARTIAL HANDS
  - BODY BUILD-UPS and RESTORATION
  - LEG & FOOT COVERS
    - ATROPHIED LEG BUILD-UPS
    - BOW LEG PADS
    - NOSE, EAR, EYE and FACIAL RESTORATIONS
    - PARTIAL FEET
    - COSMETIC FINGERS & ARMS
    - EYEBROWS
    - BREAST RESTORATIONS
    - MASKING and STRETCH HOSE
    - SUPPLIES

Realastic Industries

3010 Lyons Road

Austin, Texas 78702



READY TO FIT HAND, WRIST AND FINGER BRACES,  
CERVICAL BRACES, TRAINING AIDS AND SPECIAL  
APPLIANCES – PLUS A COMPLETE LINE OF SPINAL  
AND LEG BRACES CUSTOM MANUFACTURED TO  
YOUR MEASUREMENTS.

**COMPLETE CATALOG AVAILABLE ON REQUEST**



**C. D. DENISON** ORTHOPEDIC APPLIANCE CORP.  
220 W. 28th Street – Baltimore, Md. 21211



All Your Needs For  
ORTHOPEDIC AND PROSTHETIC APPLIANCES  
SUPPLIES

Under One Roof

*Coutils • Moleskins • Brocades • Elastics  
Non-Elastics • Nylon Lacings • Buckles  
Tools • Air Foam • Steels • Vibretta*

**Now**


**LEATHER • VELCRO®**

*New concept in fastening*


distributed to the  
Orthopedic and Prosthetic Professions  
by

***L. Laufer & Co.***


50 West 29th Street • New York 1, N.Y.




**A.** Rolling SUPPLY CART goes in elevators, wherever needed. Epoxy finish.




**B.** Orthopedic FELT & FOAM by the yd.




**C.** Velcro-held RIB BELTS, ANKLETS & WRISTLETS




**D.** LIMBO CYCLE adjusts to patient's needs, usable even with one arm or leg.




**E.** OVER BED TABLE — Chrome with tilt or non-tilt Micarta® top, adjustable ht.




**F.** Chrome BATH TUB SEAT — length adjusts.




**G.** Anodized alum. TOILET RAILS, seat or floor style.




**H.** Anodized alum. or chrome COMMODES, std. or adjustable ht.




**I.** AB TENSION PILLOW—broadcloth or vinyl.




**J.** AB INTERMITTENT TRACTION UNIT for over door or prone use.



**K.** CERVICAL COLLARS for either extension or flexion.



**L.** Adjustable width TRAPEZE FLOOR STAND



**M.** FOLDING WALK-AID, folds to 3 1/2". Anodized alum., adjustable ht. Std. models also.

Write for full descriptions, prices

**MAJOR LAB**  
MANUFACTURERS


4408 N. Sewell, Oklahoma City, Ok. 73118

# Rigidly spinal

Spinal appliances are our only product. They're sold only through ethical dispensing orthotists. So we try to make each of our collars and braces the best of its kind - like our new Two-Post Cervical Brace - for maximum acceptance by your doctors and patients. And we back them up with service to you.

Service like 24-hour delivery of prescription braces anywhere in the country. Plus a price structure that can make our service your most profitable way to fill spinal appliance prescriptions.

For more information, write P. O. Box 1299, Winter Park, Florida 32789.



**Florida Brace Corporation**

MADE TO ORDER

*Special Orthopedic Shoes*

ANY TYPE — ANY STYLE

Send us the negative casts, we make the Shoes  
to your or your Doctor's specifications.

WE SPECIALIZE IN EXTREME CASES.

For information and price list, also free measuring  
charts, write to:

**ROBERT O. PORZELT**

3457 LINCOLN AVENUE • CHICAGO 13, ILLINOIS

Telephone: BI 8-0343



### CLASSIFIED ADVERTISEMENTS

Rates: First 35 words, \$24.00 (minimum). Additional words \$1.00 each. Situations wanted advertisements half rate. Space rate additional: 13x13 picas - \$105.00; 13x18 picas - \$180.00.

Mail addressed to National Office forwarded unopened at no charge.

Classified Advertisements are to be paid in advance; checks should be made payable to "Orthotics and Prosthetics." Send to: Editor, Orthotics and Prosthetics; 1440 N St., N.W., Washington, D.C. 20005.

#### **CERTIFIED PROTHETIST-ORTHOTIST**

Will work in urban university setting doing research and training in complex experimental problems. Immediate opening, good salary, fringe benefits. Two years experience.

Mrs. Curtiss  
Box 26089  
Grady Hospital  
Atlanta, Georgia 30303

**The 1972 NATIONAL ASSEMBLY of the  
AMERICAN ORTHOTIC AND PROSTHETIC ASSOCIATION**

will be held

**OCTOBER 22-25, 1972**

at the

**STOUFFER'S RIVER FRONT  
ST. LOUIS, MISSOURI**

For Program Details and Registration Information  
write to

**The American Orthotic and Prosthetic Association  
1440 N Street, N.W.  
Washington, D. C. 20005**

The Assembly is open to all who are interested in the  
rehabilitation of the orthopedically disabled





ORTHOTICS AND PROSTHETICS  
INVITES THE SUBMISSION OF ALL ARTICLES AND MANUSCRIPTS  
WHICH CONTRIBUTE TO ORTHOTIC AND  
PROSTHETIC PRACTICE, RESEARCH, AND  
EDUCATION

All submitted manuscripts should include:

1. THE ORIGINAL MANUSCRIPT AND TWO COPIES. If possible, the duplicate manuscripts should be complete with illustrations to facilitate review and approval.
  2. BIBLIOGRAPHY. This should be arranged alphabetically and cover only references made in the body of the text.
  3. LEGENDS. List all illustration legends in order, and number to agree with illustrations.
  4. ILLUSTRATIONS. Provide any or all of the following:
    - a. Black and white glossy prints
    - b. Original drawings or charts
- Do *not* submit:
- a. Slides (colored or black & white)
  - b. Photocopies

PREPARATION OF MANUSCRIPT

1. Manuscripts must be TYPEWRITTEN, DOUBLE-SPACED and have WIDE MARGINS.
2. Indicate FOOTNOTES by means of standard symbols (\*).
3. Indicate BIBLIOGRAPHICAL REFERENCES by means of Arabic numerals in parentheses (6).
4. Write out numbers less than ten.
5. Do not number subheadings.
6. Use the word "Figure" abbreviated to indicate references to illustrations in the text (... as shown in Fig. 14)

PREPARATION OF ILLUSTRATIONS

1. Number all illustrations.
2. On the back indicate the top of each photo or chart.
3. Write author's name on the back of each illustration.
4. Do not mount prints except with rubber cement.
5. Use care with paper clips; indentations can create marks.
6. Do not write on prints; indicate number, letters, or captions on an overlay.
7. If the illustration has been published previously, provide a credit line and indicate reprint permission granted.

NOTES:

- Manuscripts are accepted for exclusive publication in ORTHOTICS AND PROSTHETICS.
- Articles and illustrations accepted for publication become the property of ORTHOTICS AND PROSTHETICS.
- Rejected manuscripts will be returned within 60 days.
- Publication of articles does not constitute endorsement of opinions and techniques.
- All materials published are copyrighted by the American Orthotic and Prosthetic Association.
- Permission to reprint is usually granted provided that appropriate credits are given.

# BRACES CAN'T BREAK DOWN THE SYSTEM.

The Alden system builds shoes that fit the foot. Accommodate prescriptions. And have the strength to accept internal and external braces.

Leather & Steel provide the strength. Leather heel bases. Heavy leather insoles. Long and strong inside counters. And tempered steel shanks constructed to take the pressure.

For complete information, write: Alden Shoe Company, Middleborough, Mass. 02346

## THE ALDEN FOOT BALANCE SYSTEM.



ALDEN SHOE COMPANY



**AMERICAN ORTHOTIC AND PROSTHETIC ASSOCIATION**  
**1440 N STREET, N.W.**  
**WASHINGTON, D.C. 20005**